

# HTM6002 Pendant

## TECHNICAL INFORMATION

Medicare pendant units are lightweight portable call devices that can be worn around the neck. The soft fabric lanyard features a strain release clip for safety.

The pendant is designed to provide an easy-to-operate tactile feel, and has a luminescent label for easy location in the dark. The pendant is manufactured to withstand the rigours of even the toughest care facility.

0800 849 5123  
www.medicare systems.co.uk

medicare  
SYSTEMS



### Enclosure

ABS moulded

### Label

Lexan with Microban® additive  
Luminescent properties for easy location in the dark

Guaranteed no fade

Can be cleaned safely with anti-bacterial cleaners

Compliant with FDA and USP requirements

### Lanyard

Nylon fabric, safety snap clip

### Programmable call type

Call

Assistance

Emergency

Low Battery

### Power

2 x AAA Alkaline batteries

### Dimensions (mm)

W x D x H: 42 x 28 x 92

### Weight including batteries

26g

### Battery life

Normal operating conditions up to 2 years

### Radio

Conforms to European EN 301 489-3 and ETSI EN 300

220-3 (for ERP @ 10mW)