



Single Button Pear Push Lead

HTM6017

The HTM6017 pear push lead is manufactured with an Anti-Microbial additive in both the enclosure and lead with built in suspension and safety snap this design is both compliant and durable.

Comes complete with replaceable linen clip for ease of use.

Operation of the device activates a normal call on the Medicare nurse call system and provides a flashing reassurance light as feedback to the user.

The device also activates an accessory call on the nurse call system when unplugged from the call point to alert staff.

Enclosure	ABS moulded manufactured with an anti-microbial additive.
Dimensions (H x W x D)	Pear Push Hand Set - 134mm x 47mm x 17mm Cable length - 2m expanding to 3m
Weight	190g
Label	Lexan with Microban additive, Luminescent for easy location in the dark, Compliant with FDA and USP requirements.
Call Facilities	Call and Lead Tamper.