



Double Button Pear Push Lead

HTM6018

The HTM6018 pear push lead is manufactured with an Anti-Microbial additive in both the enclosure and lead with built in suspension and durable design.

Operation of the device activates a normal call on the Medicare nurse call system and provides a flashing reassurance light as feedback to the user.

The device also activates an accessory call on the nurse call system when unplugged from the call point to alert staff.

When paired with a HTM6033 radio mains switch the light switch can operate either a reading or room light remotely.

Enclosure	ABS moulded manufactured with an anti-microbial additive.
Dimensions (H x W x D)	Pear Push Hand Set - 134mm x 47mm x 17mm Cable length - 2m expanding to 3m
Weight	195g
Label	Lexan with Microban additive, Luminescent for easy location in the dark, Compliant with FDA and USP requirements.
Call Facilities	Call, Lead Tamper and light switch for pairing with a HTM6033 radio mains switch.