



Infrared iButton Call Point

HTM6500iBiR

The HTM6500iBiR model is supplied with a pull cord and two sockets as standard. This model also encompasses the iButton, offering staff attendance and identification. Socket 1 is for use with a wide variety of assistive technology, such as floor mats, PIR triggers and much more. Socket 2 provides a connection to a pear push lead which comes complete with a clip and reassurance light.

Enclosure	ABS moulded c/w wall mounting bracket manufactured with an Anti-Microbial additive.
Dimensions (H x W x D)	170mm x 32mm x 96mm
Weight	251g
Label	Lexan with Microban additive, Luminescent for easy location in the dark, Compliant with FDA and USP requirements.
Radio Specifications	433.92MHz
Conformity	Conforms to European EN 301 489-3 and ETSI EN 300 220-3 (for ERP@10mW) RED
Power	2 x AA alkaline batteries
Call Facilities	Call, Assistance, Attendance, Emergency, Low battery, Accessory and Lead tamper.