



Combined Radio & Infrared Pendant

HTM6541

This Pendant is a portable nurse call product that can be used to trigger an event on the nurse call system depending on the situation of use. The soft fabric lanyard features a strain release clip for safety.

Enclosure	ABS moulded and high impact styrene.
Dimensions (H x W x D)	42mm x 28mm x 92mm
Weight	68g
Label	Lexan with Microban additive, luminescent for easy location in the dark.
Radio Specifications	433.92MHz
Conformity	ETSI EN 300 220-3 (for ERP@10mW), ETSI EN 301 489-3
Power	2 x AAA alkaline batteries.
Call Facilities	Call, Assistance, Emergency, Attack and Low battery,