



# Combined Radio & Infrared Mini Pendant

## HTM6542

This lightweight pendant has the option of being worn either around the neck on a lanyard or with a MED103 watch strap.

It is a portable nurse call product which can be programmed to transmit either a call, assistance or emergency type of event, depending on the situation of use. The soft fabric lanyard features a strain release clip for safety.

|                        |   |
|------------------------|---|
| Enclosure              | ABS moulded and manufactured with an Anti-Microbial additive.   |
| Dimensions (H x W x D) | 46mm x 63mm x 16mm  |
| Weight                 | 24g   |
| Label                  | Lexan with Microban additive, Luminescent for easy location in the dark, Compliant with FDA and USP requirements. |
| Radio Specifications   | 433.92MHz   |
| Conformity             | Conforms to European EN 301 489-3 and ETSI EN 300 220-3 (for ERP@10mW) RED  |
| Power                  | Internal rechargeable coin cell battery LIR2450   |
| Call Facilities        | Call, Assistance, Attack, Emergency and Low battery,  |